

Thermo Plate

MATS-U4020WF

The Clear Glass Heating Stage, Thermoplate MATS-U4020WF is now released for universal stereomicroscope applications.

MATS-U4020WF

... allows wide working area that is clear and heated.
... is applicable to most type of stereomicroscopes.



The glass heating stage height can be adjusted by the screws installed at the bottom of each leg.
(The glass heating stage is applicable to microscope stage height of 75mm to 100mm.)

For OLYMPUS



For NIKON



Specification

Preset Temperature

Ambient to 50 Celsius

Accuracy of Temperature Indication

+/- 0.3 Celsius

Dimension of Stage

W430 x D205 x H75-100(mm)

Dimension of Heating Glass

*Height from table top to under surface of glass heating stage frame.

Glass Thickness

W400 x D175 (mm)

1.6 (mm)

TOKAI HIT Co., Ltd.

306-1 Gendoji-cho, Fujinomiya-shi, Shizuoka-ken, Japan 418-0074

TELEPHONE: (81) 544 24 6699 FAX: (81) 544 24 6641

<http://www.tokaihit.com> E-mail: thermoplate@tokaihit.com

PL4020WFENT_20030906